

MAJOR LACERATION TRAY

EMERGENCY PROCEDURE TRAYS



CASE/12 96-1710 WITH ECONO FLOOR GRADE INSTRUMENTS

INSTRUMENTS

Curved Mosquito Hemostat
Adson Dressing Forceps
Adson Tissue Forceps
Webster Needle Holder
Straight Iris Scissors

NEEDLES/SYRINGES

Needle, 30ga x 1"
Needle, 27ga x 1 1/4"
Needle, 25ga x 5/8"
Needle, 18ga x 1 1/2"
10cc Control Syringe, L/L

ABSORBENT MATERIAL

(10) Gauze Sponges, 4 x 4, 12-ply
(5) Gauze Sponges, 2 x 2, 8-ply
(2) Blue Linen OR Towels
Fenestrated OR Towel

CONTAINERS

(2) Medicine Cups

INJECTABLES

Lidocaine, 1%, 30ml Vial

MISCELLANEOUS

Tray

MINOR LACERATION TRAY — A



CASE/12 96-1715 WITH ECONO FLOOR GRADE INSTRUMENTS

INSTRUMENTS

Webster Needle Holder
Straight Iris Scissors
Curved Mosquito Hemostat
Adson Dressing Forceps

NEEDLES/SYRINGES

Needle, 27ga x 1 1/4"
Needle, 25ga x 5/8"
Needle, 18ga x 1 1/2"
10cc Syringe, L/L

ABSORBENT MATERIAL

(10) Gauze Sponges, 4 x 4, 12-ply
(5) Gauze Sponges, 2 x 2, 12-ply
(2) Blue Linen OR Towels
Fenestrated OR Towel

CONTAINERS

(2) Medicine Cups, 2 oz.

MISCELLANEOUS

Tray

MINOR LACERATION TRAY — B



CASE/12 96-1718 WITH ECONO FLOOR GRADE INSTRUMENTS

INSTRUMENTS

OR Scissors, S/S
Baumgartner Needle Holder
Tissue Forceps

ABSORBENT MATERIAL

(10) Gauze Sponges, 4 x 4, 12-ply
(5) Gauze Sponges, 2 x 2, 12-ply
Absorbent Towel, White, 17 x 19
Fenestrated Drape, 18 x 26

MISCELLANEOUS

Tray

MINOR LACERATION TRAY



CASE/12 96-1716 WITH PLASTIC/WIRE-FORM INSTRUMENTS

INSTRUMENTS

Wire Formed Scissors
Plastic Dressing Forceps
Plastic Curved Hemostat

NEEDLES/SYRINGES

Needle, 27ga x 1 1/4"
Needle, 25ga x 5/8"
Needle, 18ga x 1 1/2"
10cc Syringe, L/L

CONTAINERS

(2) Medicine Cups, 2 oz.
Tray

ABSORBENT MATERIAL

(10) Gauze Sponges, 4 x 4, 12-ply
(5) Gauze Sponges, 2 x 2, 12-ply
Absorbent Towel, White, 17 x 19
Fenestrated Drape, 18 x 26

MINOR LACERATION TRAY

CASE/12 96-1717 WITHOUT INSTRUMENTS

FACIAL PLASTIC LACERATION TRAY

CASE/10 96-1707 WITH ECONO FLOOR GRADE INSTRUMENTS



INSTRUMENTS

Curved Mosquito Hemostat
Adson Dressing Forceps
Adson Tissue Forceps
Webster Needle Holder
Iris Scissors

CONTAINERS

(2) Medicine Cups, 2 oz.

NEEDLES/SYRINGES

Needle, 30ga x 1"
Needle, 27ga x 1 1/4"
Needle, 25ga x 5/8"
Needle, 18ga x 1 1/2"
10cc Control Syringe, L/L

MISCELLANEOUS

Tray

ABSORBENT MATERIAL

(10) Gauze Sponges, 4 x 4, 12-ply
(5) Gauze Sponges, 2 x 2, 12-ply
Fenestrated OR Towel
Blue Linen OR Towels

WOUND CLOSURE TRAY

CASE/25 96-1730 WITH ECONO FLOOR GRADE INSTRUMENTS



INSTRUMENTS

Webster Needle Holder
OR Scissors, S/B
Tissue Forceps
Curved Hemostat

MISCELLANEOUS

Tray
Refuse Bag w/ Tie

PREP

PVP Swabstick
Alcohol Swabstick

ABSORBENT MATERIAL

(4) Gauze Sponges, 4 x 4, 12-ply
Absorbent Towel, White, 17 x 19
Fenestrated Drape, 18 x 26

WOUND CLOSURE TRAY

CASE/25 96-1731 WITH PLASTIC/WIRE-FORM INSTRUMENTS



INSTRUMENTS

Wire Formed Scissors
Plastic Dressing Forceps
Plastic Curved Hemostat

MISCELLANEOUS

Tray
Refuse Bag w/ Tray

PREP

PVP Swabstick
Alcohol Swabstick

ABSORBENT MATERIAL

(4) Gauze Sponges, 4 x 4, 12-ply
Absorbent Towel, White, 17 x 19
Fenestrated Drape, 18 x 26

WOUND CLOSURE TRAY

CASE/25 96-1759 WITHOUT INSTRUMENTS

INCISION AND DRAINAGE TRAY

CASE/25 96-1713 WITH ECONO FLOOR GRADE INSTRUMENTS



INSTRUMENTS

OR Scissors, S/B
Dressing Forceps
Straight Hemostat
Scalpel, #11 Long

MISCELLANEOUS

Refuse Bag w/ Tie
Tray

PREP

PVP Swabsticks, 3's
Alcohol Swabsticks, 3's

ABSORBENT MATERIAL

(4) Gauze Sponges, 4 x 4, 12-ply
Fenestrated Drape, 18 x 26
Towel/Drape, 17 x 19, 4-ply

INCISION AND DRAINAGE TRAY

CASE/25 96-1753 WITH PLASTIC/WIRE-FORM INSTRUMENTS



INSTRUMENTS

Wire Formed Scissors
 Plastic Dressing Forceps
 Straight Plastic Hemostat
 Scalpel, #11 Long

MISCELLANEOUS

Refuse Bag w/ Tie
 Tray

PREP

PVP Swabsticks, 3's
 Alcohol Swabsticks, 3's

ABSORBENT MATERIAL

(4) Gauze Sponges, 4 x 4, 12-ply
 Fenestrated Drape, 18 x 26
 Towel/Drape, 17 x 19, 4-ply

INCISION AND DRAINAGE TRAY

CASE/25 96-1755 WITHOUT INSTRUMENTS

EPISTAXIS TRAY

CASE/25 96-1743 WITH ECONO FLOOR GRADE INSTRUMENTS



INSTRUMENTS

Plastic Nasal Speculum
 Bayonet Forceps

ABSORBENT MATERIAL

(2) Gauze Sponges, 4 x 4, 12-ply
 Nasal Packing

MISCELLANEOUS

Tray

NAIL TRAY

CASE/12 96-1838



INSTRUMENTS

Mini Scalpel
 Plastic Tissue Forceps

PREP

(2) Alcohol Prep Pads

ABSORBENT MATERIAL

(3) Gauze Sponges, 4 x 4, 12-ply
 Absorbent Towel

MISCELLANEOUS

Tray

MINOR DRESSING TRAY

CASE/25 96-1702 WITH ECONO FLOOR GRADE INSTRUMENTS



INSTRUMENTS

OR Scissors, S/B
 Dressing Forceps
 Adson Tissue Forceps

PREP

PVP Pad
 Alcohol Prep Pad

ABSORBENT MATERIAL

(3) Gauze Sponges, 4 x 4, 12-ply
 Towel/Drape, 17 x 19, 4-ply

MISCELLANEOUS

Microport Tape, 1 x 36
 Refuse Bag w/ Tie
 Tray

STANDARD PATIENT CARE TRAYS