



Small Bore

PRODUCT CODE/UPC	DESCRIPTION	FEATURES	UNITS PER CASE
MX2311L 50351688508339	3-Way Stopcock with Swivel Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50




MX2341L 50351688508322	4-Way Stopcock with Swivel Male Luer Lock, Non-Vented Blue Cap on Middle Port	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50
---------------------------	-------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----




MX5061L 50351688506847	1-Way Stopcock with Swivel Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50
---------------------------	-------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----




MX5311L 50351688506793	3-Way Stopcock with Swivel Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50
---------------------------	-------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----



MX5341L 50351688506786	4-Way Stopcock with Swivel Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50
---------------------------	-------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----




MX5341LN 50351688506779	4-Way Stopcock, Swivel Male Luer Lock, Non-Vented Red Cap on Middle Port	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50
----------------------------	--------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----




Small Bore

PRODUCT CODE/UPC	DESCRIPTION	FEATURES	UNITS PER CASE
MX5342L 50351688506762	2-Gang, 4-Way Stopcock with Swivel Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	25




Large Bore - Hi Flo™


PRODUCT CODE/UPC	DESCRIPTION	FEATURES	UNITS PER CASE
MX4312L 50351688507752	2-Gang, 3-Way Hi-Flo™ Stopcock with Extended Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	25




MX4313L 50351688507721	3-Gang, 3-Way Hi-Flo™ Stopcock with Extended Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	25
---------------------------	-------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----



MX4311L 50351688507769	3-Way Hi-Flo™ Stopcock with Extended Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50
---------------------------	-----------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----




MX4341L 50351688507646	4-Way Hi-Flo™ Stopcock with Extended Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50
---------------------------	-----------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----




Large Bore - Hi Flo™


PRODUCT CODE/UPC	DESCRIPTION	FEATURES	UNITS PER CASE
MX434B1L 50351688507585	4-Way Hi-Flo™ Stopcock with Blue Handle, Extended Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	50




MX9311L 50351688504188	3-Way Hi-Flo™ Stopcock with Swivel Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC • Yellow Handle 	50
---------------------------	---------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------	----



MX9341L 50351688504171	4-Way Hi-Flo™ Stopcock with Swivel Male Luer Lock, Non-Vented Blue Cap on Middle Port	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC • Yellow Handle 	50
---------------------------	---------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------	----



MX9342L 50351688504164	2-Gang, 4-Way Hi-Flo™ Stopcock with Swivel Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC • Yellow Handles 	25
---------------------------	-----------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----



MX4342LM 50351688507639	2-Gang, 4-Way Hi-Flo™ Stopcock Manifold, Extended Male Luer Lock	<ul style="list-style-type: none"> • Non-DEHP Formulation • Latex-Free • Lipid Resistant • Non-PVC 	25
----------------------------	------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------------------------------------------	----

