

Angiographic Accessory Products

ANGIOPLASTY ACCESSORIES

Hemostatic Y-Adapters

- Accepts a 9 Fr catheter
- Smooth turning Tuohy-Borst valve with rotating adapter
- New Myshell Lite Y-adapter is push button activated
- Available as a PTCA Accessory Kit with a guidewire insertion tool
- All items packed sterile
- Latex free

Product No.	Description	Quantity
77-400260	Hemostatic Y-Adapter - 9.5 French	25/Cs
77-400880	Hemostatic Y-Adapter - 9.0 French	25/Cs
77-400261	Hemostatic Tri-Adapter	25/Cs
77-400266	PTCA Accessory Kit (Contains: Y-Adapter & Metal Insertion Tool)	25/Cs
77-400262	PTCA Accessory Kit (Contains: Y-Adapter, Metal Insertion Tool & Torque Device)	25/Cs
77-400881	PTCA Accessory Kit (Contains: 9.0 French Y-Adapter, Metal Insertion Tool & Wide-Grip Torque Device)	25/Cs
77-400586	PTCA Accessory Kit w/S.C. (Contains: 77-400262 and 3-Way H.P. Stockcock)	25/Cs
77-400626	PTCA Accessory Kit (Contains: Y-Adapter, Metal Insertion Tool & Alligator-Yellow)	25/Cs
77-400263	Guidewire Insertion Tool, Metal	25/Cs
77-400264	Torque Device for Wires .014" - .022"	25/Cs
77-400879	Wide-Grip Torque Device .014" - .022"	25/Cs
0203NA	Alligator Torque Device .010" - .020" (Yellow)	48/Cs
77-400687	Alligator Torque Device .010" - .020" (Yellow)	10/Cs
0203NC	Alligator Torque Device .025" - .040" (Orange)	48/Cs
77-400688	Alligator Torque Device .025" - .040" (Orange)	10/Cs
0210QA	Myshell Lite Y-Connector	25/Cs
77-400590	PTCA Accessory Kit w/Myshell Lite Y-Connector (Contains: Myshell Lite, Metal Insertion Tool & Torque Device)	25/Cs



Getting to the Heart of Quality

ANGIOPLASTY ACCESSORIES



Elite™ Angioplasty Inflation Device

- Compact ergonomic design for optimal handling
- Suitable for right-and-left-handed users
- 30 ATM pressure glow gauge is reliable for all procedures
- Adjustable activation switch swivels 90° for side or top positioning and maximum comfort
- Clear 20 cc barrel is easy to fill and debubble
- Courser threading on handle for faster balloon inflation
- Latex free

Product No.	Description	Quantity
0218TK	Elite Inflation Device & 3-Way Stopcock	12/Cs
77-400843	Elite Inflation Device w/S.C. STD PTCA Kit	10/Cs
77-400842	Elite Inflation Device w/S.C. Alligator® PTCA Kit	10/Cs
77-400844	Elite Inflation Device w/S.C. Myshell Lite PTCA Kit	10/Cs