

# Angiographic Accessory Products

## CONTRAST PRODUCTS

### Contrast Saver Sets

- Minimizes contrast waste
- Maximizes potential savings
- Two-part system with or without floating disk/ball
- All items packed sterile
- Latex free

Product No.	Description	Quantity
77-400033	Contrast Saver Burette Set w/Disk	10/Cs
77-400871	60" Contrast Saver Set (POD Chamber w/Ball)	25/Cs
77-420530	72" Contrast Management Device (Cone Chamber w/Ball)	25/Cs
77-400368	In-Line Contrast Saver Set w/o Burette	25/Cs
77-400278	Mini Spike w/Vented 3" Tube & Stopcock	10/Cs
77-400065	Economy Mini Spike w/Stopcock	10/Cs
77-300733	Vented Universal Dust Caps	25/Cs
77-400208	Non-Vented Universal Dust Caps	25/Cs



### Contrast Administration Sets

- Allows smooth bubble-free priming of closed system
- Available in standard I.V. and large bore tubing
- All items packed sterile
- Latex free

Product No.	Description	Quantity
77-530042	72" Large Bore Vented Tubing	25/Cs
77-530006	72" Standard Bore Vented Tubing	25/Cs
77-530003	72" Vented Bifurcated Standard Bore	25/Cs
77-530002	72" Vented Pressurized Spike	25/Cs



### CT Transfer Set w/Swabbable Valve

- Valve is opened only when a syringe tip is inserted
- Large wings provide added grip with attaching to syringe
- Latex free

Product No.	Description	Quantity
77-3200	20" Fluid Transfer Set LL w/Wings	100/Cs
77-3208	20" Fluid Transfer Set LL w/Ball Valve on Spike	100/Cs
77-3041	6" Fluid Transfer Set w/o Swabbable Valve	100/Cs



# Getting to the Heart of Quality

## CONTRAST PRODUCTS

### Low Pressure Coiled Tubing Sets

- Single and Y-Split versions available
- Ideal for CT and MRI procedures
- Latex free

Product No.	Description	Quantity
77-400160	60" Coiled Low Pressure Tubing	50/Cs
C405-2601	60" Coiled Y-Set w/Single Check Valve	100/Cs
C405-1559	60" Straight Y-Set w/Single Check Valve	100/Cs



## FLUID ADMINISTRATION

### Fluid Administration Sets

- Allows smooth, bubble-free priming of closed system
- Available vented and non-vented
- All items packed sterile
- Latex free

Product No.	Description	Quantity
77-530000	72" Vented Standard Tubing w/Macro Drip Chamber	25/Cs
77-530001	72" Non-Vented Standard Tubing w/Macro Drip Chamber	25/Cs
77-301042	Standard Non-Vented Tubing w/.15 Micron Filter & Macro Drip Chamber	25/Cs
77-530045	72" Vented Fluid Administration Line w/Dual Check Value	25/Cs
77-530007	72" Non-Vented Bifurcated w/Macro Drip Chamber	25/Cs
77-400834	108" Universal Admin Set w/Macro Drip Chamber	25/Cs
77-400837	72" Red Stripe Non-Vented w/Macro Drip Chamber	25/Cs
77-400836	72" Blue Stripe Non-Vented w/Macro Drip Chamber	25/Cs
77-400835	72" Green Stripe Non-Vented w/Macro Drip Chamber	25/Cs
77-400838	72" Yellow Stripe Non-Vented w/Macro Drip Chamber	25/Cs



### Fluid Transfer Decanter

- Ideal for bags of saline
- Latex free

Product No.	Description	Quantity
77-400274	Decanter	50/Cs



# Getting to the Heart of Quality

## SYRINGES



### Power Injector Syringes

- CT Tri-Pack fits most injectors with free adapter
- Power Injector Syringes for most popular injectors
- Value priced
- All items packed sterile
- Latex free

Product No.	Description	Fits These Machines	Quantity
77-150193	150 ml CT Tri-Pack w/8" Fill Tube & 60" Coiled Line	Leibel-Flarsheim Angiomat® 6000, Illumena®**	50/Cs
77-301380	200 ml CT Tri-Pack w/8" Fill Tube & 60" Coiled Line (Turret Load)	Medrad MCT™, MCT Plus™, Medrad 200 Series	50/Cs
77-200590	200 ml CT Tri-Pack w/8" Fill Tube & 60" Coiled Line (Front Load)	Medrad MCT™*, MCT Plus™ Front Load*, EnVision CT™ Front Load*, Angiomat® CT9000* & CT9000 ADV*	50/Cs
77-400040	100 cc Power Injector Syringe w/J-Straw	Leibel-Flarsheim Angiomat® 3000 Injectors	10/Cs
77-400130	130 cc Power Injector Syringe	Leibel-Flarsheim Angiomat® 6000, Medrad Mark II™, II & IV, Medrad CT™ 100 Series	50/Cs
77-400230	130 cc Power Injector Syringe w/J-Straw	Leibel-Flarsheim Angiomat® 6000, Medrad Mark II™, II & IV, Medrad CT™ 100 Series	50/Cs
77-400277	150 cc Power Injector Syringe	Medrad Mark V™, Mark V Plus™, Provis®	50/Cs
77-400250	150 cc Power Injector Syringe w/J-Straw	Medrad Mark V™, Mark V Plus™, Provis®	50/Cs
77-150269	150 cc Power Injector Syringe w/J-Straw	Leibel-Flarsheim Angiomat® 6000, Illumena®**	50/Cs
77-400220	200 cc Power Injector syringe w/J-Straw	Medrad Mark V™, Mark V Plus™, Provis®	50/Cs
77-400041	8" Fill Straw		50/Cs
77-400160	60" Female to Male Low Pressure Coiled Line		50/Cs

\* Adapter required and provided free of charge. \*\* Pressure jack required and provided free of charge.

\*\*\* Angiomat and Illumena are registered trademarks of Leibel-Flarsheim; MCT, MCT Plus, EnVision CT, Mark IV, Mark V, Mark V Plus, and Provis are trademarks of Medrad.

# Angiographic Accessory Products

## SYRINGES

### Control Syringes

- Smooth operation
- Clear, glass-like polycarbonate for complete clarity
- New 8 ml and 10 ml narrow barrel syringes that are ideal for use with smaller catheters
- Available with .5 ml reservoir on plunger
- All items packed sterile
- Latex free

#### Narrow Barrel

Product No.	Description	Quantity
77-400391	8 ml Male Luer Lock Thumb Ring & Finger Rings	25/Cs
77-400392	8 ml Male Rotator Thumb Ring & Finger Rings	25/Cs
77-712391	10 ml Male Luer Lock Thumb Ring & Finger Rings	25/Cs

#### Standard Barrel

Product No.	Description	Quantity
77-300741	10 ml Male Luer Lock Thumb Ring & Finger Rings	25/Cs
77-300744	10 ml Male Rotator Thumb Ring & Finger Rings	25/Cs
77-300743	10 ml Male Luer Lock Palm Pad & Wings	25/Cs
77-300742	10 ml Male Rotator Palm Pad & Wings	25/Cs
77-301291	12 ml Male Luer Lock Thumb Ring & Finger Rings	25/Cs
77-301294	12 ml Male Rotator Thumb Ring & Finger Rings	25/Cs
77-301293	12 ml Male Luer Lock Palm Pad & Wings	25/Cs
77-301292	12 ml Male Rotator Palm Pad & Wings	25/Cs



### Colored Piston Specialty Syringes

- Smooth operation
- Ideal for fluid recognition
- Add labels to customize as needed
- All items packed sterile
- Latex free

Product No.	Description	Quantity
77-400100	10 cc Red Specialty Syringe	50/Cs
77-400101	10 cc Blue Specialty Syringe	50/Cs
77-400102	10 cc Yellow Specialty Syringe	50/Cs
77-400103	10 cc Green Specialty Syringe	50/Cs
77-400104	10 cc White Specialty Syringe	50/Cs
77-400105	20 cc Blue Specialty Syringe	50/Cs
77-400106	20 cc Green Specialty Syringe	50/Cs
77-400107	20 cc Red Specialty Syringe	50/Cs
77-400108	20 cc White Specialty Syringe	50/Cs
77-400109	20 cc Yellow Specialty Syringe	50/Cs

