

Getting to the Heart of Quality

ANGIOGRAPHY ESSENTIALS

Diagnostic Guidewires

- Straighten with finger pressure
- Produced by leading manufacturers for DeRoyal
- Packaged sterile in dispenser
- Latex free

Product No.	Description	Quantity
77-8546	J3 Fixed-Core 150 cm - .035 GUID, PTFE, FS	25/Cs
77-400485	J3 Fixed-Core 150 cm - .038 GUID, PTFE, FS	25/Cs



Angiography Percutaneous Guidewire Introducer Needles

- Sharp cannulas for low penetration force
- Latex free

Product No.	Description	Quantity
77-11606	18 Gauge 2¾" One Part w/o Wings	50/Cs
77-11078	18 Gauge 2¾" One Part w/Wings	50/Cs
77-400507	18 Gauge 4" Two Part w/Wings	25/Cs




Sharps Safety Guidewire Introducer Needle

- Assists clinicians in improving safety
- Innovative needle safety handle that advances from hub to cover needle point
- Sheath that covers needle shaft simultaneously advances as the handle is activated
- Provides protection from blood on shaft as well
- Promotes compliance with H.R. 5178, the Needlestick Safety and Prevention Act
- Meets OSHA BBP guidelines
- Latex free

Product No.	Description	Quantity
SPG-101-18-07	18 Gauge x 7 cm SPpro™ GWI Sharps Safety Needle	25/Cs
SPG-101-19-07	19 Gauge x 7 cm SPpro™ GWI Sharps Safety Needle	25/Cs
SPG-101-21-07	21 Gauge x 7 cm SPpro™ GWI Sharps Safety Needle	25/Cs

*SPpro is a trademark of NeoMed.



Angiographic Accessory Products

WASTE COLLECTION

Waste Collection Stations

- Facilitates the safe collection of fluid waste
- Also allows the administration of saline and contrast solutions
- Minimizes contact with infectious blood and bodily fluids
- Leak resistant dual check valves make disposal clean and easy
- A wide variety of waste collection stations
- Latex free unless otherwise noted
- All items packed sterile

Product No.	Description	Quantity
77-300695	Waste Station w/Macro Drip Chamber, Dual Check Valve & Drape Clip	25/Cs
77-300692	Waste Station w/Vented Spike, Dual Check Valve & Drape Clip	25/Cs
77-300691	Waste Station w/Vented Spike, Dual Check Valve, Anti-Reflux & Drape Clip	25/Cs
77-400650	Waste Station w/Drip Chamber Plus Contrast Line	20/Cs
77-400215	Waste Station Plus w/Blue Stripe Contrast Line, Dual Check Valve & Drape Clip	20/Cs
77-300735	Waste Station Plus Blue Contrast Line & Drape Clip	20/Cs
77-400192	Waste Station Plus In-Line Contrast Saver (Use w/Mini Spike, Pg. 11)	10/Cs
77-300726	Abscess Evacuation Set	25/Cs
77-300752	Waste Bag w/60" Tubing & Male Luer Lock	25/Cs
77-300783	Waste Bag w/60" Tubing, In-Line Check Valve & Male Luer Lock	25/Cs



SplashStop™

- A safe alternative to messy, open basins
- Closed bowl minimizes staff and patient exposure to blood-borne pathogens
- A true engineering and work practice control
- Secure, non-splatter lid with tapered insertion site
- Safe, one-step biohazard warning label
- All items packed sterile
- Latex free

Product No.	Description	Quantity
77-301370	SplashStop Waste Collection Bowl	20/Cs
77-301371	SplashStop Waste Collection Bowl w/Solidifier	20/Cs

* SplashStop is a trademark of DeRoyal.

