

Angiographic Accessory Products

ANGIOGRAPHY ESSENTIALS

Angiographic Equipment Covers

- Clear polyethylene
- Ensures sterility and protects equipment during procedures
- Many sizes and styles available in addition to those listed*
- All items packed sterile
- Latex free

Product No.	Description	Quantity
77-100011	Round 14" Diameter	25/Cs
77-100012	Round 22" Diameter	25/Cs
77-100013	Round 27" Diameter	25/Cs
77-103030T	30" x 30" w/Tape	25/Cs
77-104030T	40" x 30" w/Tape	25/Cs
77-104060T	40" x 60" w/Tape	25/Cs

*Contact your DeRoyal representative with the dimensions you require.



Guidewire Bowl w/Lid

- Patented Design with built up lip and recessed bottom
- Designed to hold Guidewire in place and allow easy removal
- All items packed sterile
- Latex free

Product No.	Description	Quantity
32-1345	Guidewire Bowl w/Lid	25/Cs



Rotating Adapters

- Allows smooth operation for catheter adjustments to prevent torque build-up
- Made of high-quality polycarbonate
- All items packed sterile
- Latex free

Product No.	Description	Quantity
77-400258	Female to Male Rotating Adapter	25/Cs
77-400259	Male to Male Rotating Adapter	25/Cs



Getting to the Heart of Quality



ANGIOGRAPHY ESSENTIALS

Labels For Angiography & Cath Labs

- Meet 2006 Joint Commission (JCHAO) requirements for labeling:
 - All medications and solutions
 - On and off sterile field
 - When medication or solution is transferred from original packaging to another container
 - In all hospitals, critical access facilities, ambulatory care settings and office based-surgery
- Standard kits with 1/2" x 1-3/4" labels, color-coded or white
- Clear, bold printing in black letters
- Most commonly used angiography medications/solutions plus blank labels
- With or without utility markers
- Custom label kits available to fit individual procedures or facility requirements
- Latex free

Product No.	Description	Quantity
77-400563	Colored Label Kit* 1/2" x 1-3/4"	14/Pk, 100Pks/Cs
77-400564	Colored Label Kit* 1/2" x 1-3/4" w/Marker	14/Pk, 100Pks/Cs
77-400565	White Label Kit 1/2" x 1-3/4"	14/Pk, 100Pks/Cs
77-400566	White Label Kit 1/2" x 1-3/4" w/Marker	14/Pk, 100Pks/Cs

*Colored Label Kit contains: 4 Blue Saline, 2 Green Contrast, 2 Red Nitro, 2 Yellow Lidocaine, 2 White Heparin, 2 White Blank



Surgical Utility Markers

- Single use, utility markers
- For use on instruments, specimens, vials & tubes, basins, bowls, and other objects
- Non-smearing and fast-drying
- Fine-tip or wide-tip
- Available non-sterile or sterile for utilization in sterile field
- May be ordered separately or with kits
- Non-toxic, non-solvent, alcohol based, acid free
- Latex free

Product No.	Description	Quantity
26-0168	Surgical Utility Marker, Black, Fine Tip	50/Cs
26-0168NS	Surgical Utility Marker, Black, Non-Sterile	200/Cs
26-0168W	Surgical Utility Marker, Black, Wide Tip	50/Cs
26-0168WNS	Surgical Utility Marker, Black, Wide Tip, Non-Sterile	200/Cs
26-0168WL	Surgical Utility Marker, Black, Wide Tip, w/2 Labels 1"x3"	50/Cs

