

Angiographic Accessory Products

ANGIOGRAPHY ESSENTIALS



mBrace™ Wrist Support

- Flat front edge helps prevent rolling and provides a stable platform to support the wrist
- Soft fabric helps wick away moisture and provides comfort for the patient, especially for longer arterial line deployments
- Consistent 40° wrist extension angle from one patient to the next
- Foam pad can be used to easily adjust angle of extension and cushion IV lines on the dorsal surface of hand or wrist
- Straps adjust easily to accommodate virtually any size hand
- Compatible with radial compression devices for post-catheterization hemostasis and specially designed for use with the RadAR™ Vascular Compression Device
- Non-sterile
- Latex free

Product No.	Description	Quantity
140-0200-00	mBrace Wrist Support	30/Cs
140-0200-00S	mBrace Wrist Support, Sterile	30/Cs

* mBrace is a trademark of Semier Technologies, Inc.

Getting to the Heart of Quality

ANGIOGRAPHY ESSENTIALS



Flexible I.V. Armboard

- Malleable board for custom positioning with soft, foam cover
- Adjustable foam straps eliminate the need for skin-irritating tape
- Straps attach to board to limit arm and hand movement during IV therapy
- Can be used for extending and stabilizing the wrist during vascular access procedures, including arterial line insertions and cardiac catheterizations using the radial artery or vein
- Helps prevent patients from flexing their wrist during extended deployments for arterial lines
- Latex free

Product No.	Description	Quantity
M8125-A	Adult Flexible Arm Board (9" long x 4" wide)	25/Cs
M8125-P	Preemie Flexible Arm Board (2" long x 1" wide)	30/Cs
M8125-I	Infant Flexible Arm Board (4" long x 1" wide)	30/Cs
M-8125-C	Child Flexible Arm Board (6" long x 2" wide)	30/Cs