



Needles

One-Part Seldinger Needles Without Baseplate with Plug

Cardiac Catheterization Needles that are designed with a transparent polycarbonate hub and a very sharp "B" type arterial bevel. Our clear polycarbonate construction coupled with our smooth glide funnel allows for trouble-free guide wire insertion.



31-1815 shown



31-1827-P shown

ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		

18G NEEDLE

31-1815	Pink	18g	3.8	1-1/2	.035/.038	No
31-1822	Pink	18g	5.7	2-1/4	.035/.038	No
31-1827	Pink	18g	7.0	2-3/4	.035/.038	No
31-1835	Pink	18g	9.0	3-1/2	.035/.038	No

Tip: 

SPECIAL LENGTHS

31-1845	Clear	18g	11.0	4-1/2	.035/.038	No
31-1860	Clear	18g	15.0	6.0	.035/.038	No
31-1880	Clear	18g	20.0	8.0	.035/.038	No

Tip: 

18G NEEDLE WITH POROUS PLUG

31-1815-P	Pink	18g	3.8	1-1/2	.035/.038	No
31-1822-P	Pink	18g	5.7	2-1/4	.035/.038	No
31-1827-P	Pink	18g	7.0	2-3/4	.035/.038	No
31-1835-P	Pink	18g	9.0	3-1/2	.035/.038	No

Tip: 

Special Note: A finger guide is embedded on one side of the hub to indicate bevel placement.



Needles



One-Part Seldinger Needles Without Baseplate
(continued)

Cardiac Catheterization Needles that are designed with a transparent polycarbonate hub and a very sharp "B" type arterial bevel. Our clear polycarbonate construction coupled with our smooth glide funnel allows for trouble-free guide wire insertion.



Picture shown with porous plug

ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		

19G NEEDLE

329-355	Brown	19g	5.7	2-1/4	.032/.035	No
329-357	Brown	19g	7.0	2-3/4	.032/.035	No

Tip: 

19G NEEDLE WITH POROUS PLUG

329-355-P	Brown	19g	5.7	2-1/4	.032/.035	No
329-357-P	Brown	19g	7.0	2-3/4	.032/.035	No

Tip: 

PEDIATRIC NEEDLES: 21G NEEDLE

31-2115	Green	21g	3.8	1-1/2	.018/.021	No
---------	-------	-----	-----	-------	-----------	----

Tip: 

PEDIATRIC NEEDLES: 21G NEEDLE WITH POROUS PLUG

31-2115-P	Green	21g	3.8	1-1/2	.018/.021	No
-----------	-------	-----	-----	-------	-----------	----

Tip: 



Needles

One-Part Seldinger Needles with Courmand Plate

Cardiac Catheterization Needles are designed with a "Courmand Wing Plate" transparent polycarbonate hub and a very sharp "B" type arterial bevel. Our clear polycarbonate construction coupled with our smooth glide funnel allows for trouble-free guide wire insertion.



ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		

BSDN 18G NEEDLE

21-1815	Pink	18g	3.8	1-1/2	.035/.038	No
21-1822	Pink	18g	5.7	2-1/4	.035/.038	No
21-1827	Pink	18g	7.0	2-3/4	.035/.038	No
21-1835	Pink	18g	9.0	3-1/2	.035/.038	No

Tip: 

18G NEEDLE WITH POROUS PLUG

21-1815-P	Pink	18g	3.8	1-1/2	.035/.038	No
21-1822-P	Pink	18g	5.7	2-1/4	.035/.038	No
21-1827-P	Pink	18g	7.0	2-3/4	.035/.038	No
21-1835-P	Pink	18g	9.0	3-1/2	.035/.038	No

Tip: 

19G NEEDLE

229-355	Brown	19g	5.7	2-1/4	.032/.035	No
229-357	Brown	19g	7.0	2-3/4	.032/.035	No

Tip: 

19G NEEDLE WITH POROUS PLUG

229-355-P	Brown	19g	5.7	2-1/4	.032/.035	No
229-357-P	Brown	19g	7.0	2-3/4	.032/.035	No

Tip: 



Seldinger Needles

SELDINGER "ARTERIAL" NEEDLES – Solid Stylet

Disposable access needle typically used in the introduction of a guide wire into the vasculature. Our clear polycarbonate construction coupled with smooth glide funnel allows for trouble-free guide wire insertion.



ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
55-1829	Pink	18g	7.0	2-3/4	.035/.038	No
55-1828	Pink	18g	7.0	2-3/4	.035/.038	Yes

Tip (solid):



MODIFIED SELDINGER "ARTERIAL" NEEDLES – Hollow Stylet

Modified Seldinger Needle is a blend of the Potts Courmand point configuration and the Original Seldinger Needle large plate. For that blending reason, the Modified Seldinger are also known as Hollow Stylet Seldinger Needles. Our clear polycarbonate construction coupled with smooth glide funnel allows for trouble-free guide wire insertion.



Shown here with obturator

ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
50-1822	Pink	18g	5.7	2-5/64	.035/.038	No
50-1823	Pink	18g	5.7	2-5/64	.035/.038	Yes
50-1827	Pink	18g	7.0	2-3/4	.035/.038	No
50-1828	Pink	18g	7.0	2-3/4	.035/.038	Yes
50-1830	Pink	18g	7.6	3.0	.035/.038	No
50-1835	Pink	18g	9.0	3-1/2	.035/.038	No
50-1840	Pink	18g	10.2	4.0	.035/.038	No

Tip:



Seldinger Needles



ULTRA-MODIFIED SELDINGER NEEDLES

The Ultra Modified Seldinger Needle has an ultra-thin 19 gauge outer cannula fitted to a hollow cutting stylet. This needle is designed to make the smallest possible percutaneous puncture for delivery of a .038 guide wire. Our clear polycarbonate construction coupled with our smooth glide funnel allows for trouble-free guide wire insertion.



ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
58-0387	Brown	19g	7.0	2-3/4	.035/.038	No

Tip:



MODIFIED SELDINGER NEEDLE WITHOUT PLATE

The .035 Modified Seldinger Needle without plate has an extra-thin 19 gauge outer cannula fitted to a hollow cutting stylet. This needle is very similar to the Modified Seldinger series but without a large Kidney Base Plate. This needle will accept an .035 guide wire. The hub is clear with the stylet color coded brown denoting 19 gauge.



ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
319-0357	Brown	19g	7.0	2-3/4	.035	No

Tip:





Needles

One-Part Seldinger Needles with Kidney Plate

Cardiac Catheterization Needles are designed with a Large Kidney Plate with a transparent polycarbonate hub and a very sharp "B" type Arterial Bevel. Our clear polycarbonate construction coupled with our smooth glide funnel allows for trouble-free guide wire insertion.



18G NEEDLE

ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
51-1815	Pink	18g	3.8	1-1/2	.035/.038	No
51-1822	Pink	18g	5.7	2-1/4	.035/.038	No
51-1827	Pink	18g	7.0	2-3/4	.035/.038	No
51-1835	Pink	18g	9.0	3-1/2	.035/.038	No

Tip:



18G NEEDLE WITH POROUS PLUG

ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
51-1815-P	Pink	18g	3.8	1-1/2	.035/.038	No
51-1822-P	Pink	18g	5.7	2-1/4	.035/.038	No
51-1827-P	Pink	18g	7.0	2-3/4	.035/.038	No
51-1835-P	Pink	18g	9.0	3-1/2	.035/.038	No

Tip:



Needles



Ultra Thin Catheterization Needles One-Part

U-KSDN CARDIAC SPECIAL GAUGE

These Ultra Thin needles offer the advantage of the smallest puncture possible to deliver an .038 wire. Our clear polycarbonate construction coupled with our smooth glide funnel allows for trouble-free guide wire insertion.



Picture shown with porous plug

ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
57-0387	Brown	19g	7.0	2-3/4	.035/.038	No

Tip: 

WITH POROUS PLUG

57-0387-P	Brown	19g	7.0	2-3/4	.035/.038	No
-----------	-------	-----	-----	-------	-----------	----

Tip: 

U-SDN CARDIAC SPECIAL GAUGE

These Ultra Thin needles offer the advantage of the smallest puncture possible to deliver an .038 wire. Our clear polycarbonate construction coupled with our smooth glide funnel allows for trouble-free guide wire insertion.



Picture shown with porous plug

ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
37-0387	Clear	19g	7.0	2-3/4	.035/.038	No

Tip: 

WITH POROUS PLUG

37-0387-P	Clear	19g	7.0	2-3/4	.035/.038	No
-----------	-------	-----	-----	-------	-----------	----

Tip: 