

## Needles



### Cournand-Style Needles

#### POTTS COURNAND NEEDLES

Our Potts Cournand type needle maintains the integrity of the original design, but we have upgraded the plastic to a clear (non-shadow producing) polycarbonate. This model allows the shaft to protrude out the hub.




ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
24-1821	Pink	18g	5.7	2-5/64	.035/.038	No
23-1821	Pink	18g	5.7	2-5/64	.035/.038	Yes

Tip: 

### Cournand-Style Needles

#### MODIFIED POTTS COURNAND NEEDLES

Disposable needle set used for Carotid Angiography as well as direct arterial pressure monitoring, blood sampling and percutaneous catheterization. Our clear polycarbonate construction coupled with smooth glide funnel allows for trouble-free guidewire insertion.



Tip: 

ReOrder #	Hub Color	Gauge Size	Length		Accepts Guide	Obturator
			Cm.	In.		
20-1822	Pink	18g	5.7	2-5/64	.035/.038	No
20-1823	Pink	18g	5.7	2-5/64	.035/.038	Yes
20-1827	Pink	18g	7.0	2-3/4	.035/.038	No
20-1828	Pink	18g	7.0	2-3/4	.035/.038	Yes
20-1830	Pink	18g	7.6	3.0	.035/.038	No
20-1835	Pink	18g	9.0	3-1/2	.035/.038	No
20-1840	Pink	18g	10.2	4.0	.035/.038	No

#### 19G NEEDLES

219-0355	Brown	19g	5.7	2-5/64	.032/.035	No
219-0357	Brown	19g	7.0	2-3/4	.032/.035	No

Tip: 




---

**Needles**

---

**Chiba General Purpose Aspiration Cytology Biopsy Needle**



**CHIBA 20 GAUGE – UNMARKED**

Order	Hub Color	Stylet Color	Gauge Size	Length		Accepts Guide	Bevel Tip EchoGenic	Cannula Marked
				Cm	In			
8220	Yellow	N/A	20g	11.5	4.5	.021	No	No

Tip: 



**CHIBA 22 GAUGE - MARKED - ECHO GENIC**

Order	Hub Color	Stylet Color	Gauge Size	Length		Accepts Guide	Bevel Tip EchoGenic	Cannula Marked
				Cm	In			
8204	Clear	Black	22g	11.5	4.5	.018	Yes	Yes

Tip: 