

# DEXMED

MEDICAL PRODUCTS



X-RAY PROTECTION

**DEXMED**INC. 2008 Rt. 37 East, Suite 13, Toms River, NJ 08753  
732-831-0507 fax 732-831-1607 dexmed.com

## **X-Ray Protection Services**

### **Lead Aprons and garments**

New Lead aprons  
Custom sized Lead aprons  
Apron repair and revisions  
Apron recoveries  
Arm shielding  
Leg shielding

### **Lead apron accessories**

Mobile Lead Barriers  
Apron racks and hangers

### **Lead glasses**

Prescription Lead glasses  
Lead glasses lens replacements  
Frame replacement and repairs  
Lead / Laser protective eyewear

### **Lead gloves**

Sterile surgical gloves  
Reusable gloves and mitts  
Positioning gloves

### **Laser shielding eyewear**

Laser protective eyewear  
Lead / Laser combination eyewear



**VELCRO APRON**



**DROP AWAY APRON**



**JACKET APRON**



**TWO PIECE APRON**

## **LEAD** OPTIONS

**ULTRALITE:** This is the lightest lead you can buy unless you prefer to sacrifice protection. Our Ultralite aprons are approximately 10% lighter than our Featherweight aprons.

**FEATHERWEIGHT:** This lead is approximately 15% lighter and is more flexible than our regular lead but provides the same protection up to 100 kVp.

**REGULAR:** When economy is your goal, this option is for you. This lead has been reformulated to make it more flexible and to reduce the weight without sacrificing the protection.



CONSIDER THE  
**POSSIBILITIES**

**Fabrics, Colors Embroidery,  
Apron Styles, Custom Sizing**

Dexmed has more choices for the customer than any other manufacturer. We offer at least six types of fabrics with more than 50 colors and more than 60 choices of logos at no additional charge. Custom sizing is included, and your fit is guaranteed. In addition to the styles shown, you name it and we make it: from Economy to Designer.

**CALL DEXMED OR OUR  
LOCAL SALES REPRESENTATIVE  
FOR MORE INFORMATION**



**STRETCHYBACK APRON**



**TIE APRON**



**BUCKLE APRON**



**CUMMERBUND APRON**

### Cumberbund Apron

Back discomfort is a complaint of the past when you wear our most comfortable apron. The wide belt takes the weight off the shoulders and provides lower back support. The shoulder pads are standard for additional comfort. Should drop-away optional.

[Size chart click here](#)

[Embroidery options click here](#)



### Slot Vest Apron

Poor flexibility is a complaint of the past when you wear our most maneuverable apron. Two piece apron. The top features an adjustable stretch back with velcro shoulder and side closures. Adjustable velcro closure skirt.



**\*\* Available with Intergral Back Support**