

APRON ACCESSORIES

Thyroid Shields

Standard Thyroid

- The only thyroid made with 2" loop velcro and 1" hook, for easier closure.
- The softest, most comfortable thyroid collar made. No binding to irritate the neck.



- TS-1 Standard
- TS-2 Customer Lead
- TS-3 Ultimate
- TS-4 X-Large
- TS-5 Ring Attached
- TS-6 Permanent Attached

Ring Attached

- Keep your collar attached to your apron, or detach it with this clever design.



GS-404 (Gonad Shield Set)

This set includes 2 small, 1 medium, 1 large gonad shields and a belt all with Velcro for an interchangeable system
Gonad Shield with Belt System



WOW BELT

- A patented pulley system cinches the belt for superior lumbar support.
- The perfect addition to keep the weight off the shoulders.



SIZES

60122	25" - 30"
60148	30" - 35"
60137	35" - 40"
60155	40" - 45"
60170	45" - 50"

Cummerbund Belt

ABXS	to 25"
ABS	to 24" - 34"
ABM	to 34" - 40"
ABL	to 36" - 50"
ABX	to 40" - 60"

Lap Pad

Cat. No. LP-100

- Standard size 18" x 24"
- Protection for any area of the body.
- Any custom size available



Gonad Shields

Individual Shields: (with or without ties)

- GS-101 Small 6 x 4"
- GS-102 Medium 8 x 6"
- GS-103 Large 10 x 7"
- GS-104 Set of Four
- * Add T to Cat. No. for ties (GS-103T)

A padded support belt for extra comfort. It can be snapped, sewn on or attached by our standard loop system.



APRON ACCESSORIES

APRON ACCESSORIES

Demi-Apron Set/ Rack



- Attached to soft steel spring for ease of application and to hold securely in place.
- Aprons may be purchased separately or as a set with or without the rack.

Cat. No.	Size	Width	Length
803-S	Small	8-1/2"	10"
803-M	Medium	12"	12"
803-L	Large	16"	14"
803-X	Xlarge	20"	16"

830 4 Demi Aprons with Rack

830-R Demi Aprons Rack Only

Half Aprons

Cat. No. HA-1 to 4
(shown above)

- Velcro closure
- Buckle upon request
- Custom sizes available

		(l) x (w)
HA-1	Small	12" x 14"
HA-2	Medium	14" x 16"
HA-3	Large	16" x 18"
HA-4	Xlarge	20" x 18"

Name Tag

4 1/2" x 2" on soft Velcro.
Embroider with 12 characters.

Suspenders

When hips and waist are on the same size, you can hold the skirt up with suspenders.

Drop-Away Shoulders

Convert your apron to a Drop Away by closing the shoulders with velcro rather than sewing.

Scoliosis Shawl

Designed for use on spinal radiographic procedures involving upright or recumbent positions.



Apron Sleeves

- Attached with snaps and velcro combination for security and comfort.
- Optional permanent attachment.

SS-100	Small	9" long
SS-200	Medium	11" long
SS-300	Large	13" long
SS-500	Pair	



Elastic in waist

For additional support or a small waist, we sew 6" wide elastic in the back panel of a skirt