

SHIELDS


PORTA SHIELD Cat. No. 1001

- Very useful in helping shield areas of a patient's body during some radiographic procedures.
- Adjustable from 36" to 60" high.
- The Porta Shield is supported on a stable five legged base that moves easily with the touch of a finger.
- Covered in nylon reinforced vinyl.
- Shield size 24" x 24"



SCATTER SHIELD Cat. No. 8055 (Complete set Navy Blue) 8056 (Slide Rail Bracket) 8057 (Replacement Drapes) 8058 (L-Arm Support) 8055-C (Custom set – Call for availability)

- Protects physicians and assistants from scatter radiation from the patient and the table top. It also provides lower extremity protection for personnel in close proximity to the tube, especially during lateral head studies.
- Constructed of stainless steel and aluminum (metal parts), the Scatter Shield is perfect for use in Cardiac Cath, Special Studies, O.R., Anesthesia, EPS and Endoscopy.
- Easily mounts on any table that has a side-rail. Shield system moves with the table for uniform, continuous protection during the case.
- Foot controls are easily accessible. Physicians can easily adjust position of the Scatter Shield during a case with turnknob.





Mobile Leaded Barrier LB-2430

Overall Size: 75(h) x 30 in (w)

Window pane lead equivalent: 24 (h) x 30 in (w) 2.0mm Lead Equivalency

Opaque panel lead equivalent: 48(h) x 30 in (w) 2.0mm Lead Equivalency

Shipping weight approx: 110lb

\$3,790.00 MSRP **\$2,945.00 DEXMED**



Mobile Leaded Barrier LB-3060

Overall Size: 75(h) x 30 in (w)

Window pane lead equivalent: 60(h) x 30 in (w) 2.0mm Lead Equivalency

Opaque panel lead equivalent: 12(h) x 30 in (w) 2.0mm Lead Equivalency

Shipping weight approx: 113lb

\$5,260.00 MSRP **\$4,123.00 DEXMED**



Mobile Leaded Barrier LB-3648

Overall Size: 75(h) x 48 in (w)

Window pane lead equivalent: 36(h) x 48 in (w) 2.0mm Lead Equivalency

Opaque panel lead equivalent: 36(h) x 48 in (w) 2.0mm Lead Equivalency

Shipping weight approx: 225lb

\$6,340 MSRP **\$5,445.00 DEXMED**



Mobile Leaded Barrier LB-3672

Overall Size: 75(h) x 72 in (w)

Window pane lead equivalent: 36(h) x 72 in (w) 2.0mm Lead Equivalency

Opaque panel lead equivalent: 36(h) x 72 in (w) 2.0mm Lead Equivalency

Shipping weight approx: 320lb

\$7,820.00 MSRP **\$6,415.00 DEXMED**



Opaque Mobile Leaded Barriers LB-1012-

All Opaque Barriers have 10"x12" Windows

	MSRP	DEXMED
LB-1012-2472 2.0mm Pb, 24"x72"	\$1,980.00	\$1,782.00
LB-1012-2484 2.0mm Pb, 24"x84"	\$2,130.00	\$1,917.00
LB-1012-3072 2.0mm Pb, 30"x72"	\$2,150.00	\$1,935.00
LB-1012-3084 2.0mm Pb, 30"x84"	\$2,360.00	\$2,124.00
LB-1012-3672 2.0mm Pb, 36"x72"	\$2,300.00	\$2,070.00
LB-1012-3684 2.0mm Pb, 36"x84"	\$2,500.00	\$2,250.00
LB-1012-4572 2.0mm Pb, 45"x72"	\$2,400.00	\$2,160.00
LB-1012-4584 2.0mm Pb, 45"x84"	\$2,600.00	\$2,340.00