

DEXMED_{Inc.}

A Unique Source for All your Invasive Cardiology and Radiology Lab Needs.

Call Toll Free **1.800.8DEXMED****News****Contact / Ordering****Product line**

Sterile Packs and Kits
 Guide Wires
 Introducer Needles
 Manifold Kits
 Tubing
 Monitoring & Transducers
 Pressure Infusers
 Lead Aprons
Lead Glasses
 Laser Glasses
 Medical & Dental Surgical Loupes
 Lead Accessories
 Angiographic Syringes
 Surgical Instruments

Lead Glasses*Mouseover colors for enlargements***Escapades**

Description: The 808 is a comfortable and lightweight yet sturdy unifit nylon frame. Don't let the sharp-looking streamlined design fool you, this frame is all function. Rubber nose pads and rubberized temple bars provide comfortable grip while the wrap around design ensures that your eyes are protected. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Red; Silver; Black**Order No:** #RG-808**Biker**

Description: The 208 is a stylish wraparound frame made out of TR-90 impact resistant plastic. This wide-fit frame is perfect for larger men or anyone who has constant issues with frames being too narrow. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Silver; Red; Taupe; Matte Black; Blue/Yellow; Grey/Purple**Order No:** #RG-208**Biker Petite**

Description: Our Biker Petite is a clone of our #208 Biker in a smaller size suitable for petite and medium faces. Along with soft rubber temple pads, we've added soft rubber nose pads for the ultimate in comfort and fit. This very popular frame features our standard 0.75mm Lead equivalency Schott SF-6 HT Radiation safety lenses.

Color: Matte Black; Red; Silver**Order No:** #EL201**Vader**

Description: The Model 375 offers superior protection in a premium frame. Centered on protection, style, and performance, the frame features a head-hugging temple design with excellent splash and side protection as well as a designer flair for style. The frame is constructed out of a durable TR-90 nylon plastic for high impact and chemical resistance. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses. Excellent choice for smaller to medium faces.

Color: Black; Tortoise**Order No:** #RG-375**Maxi**

Description: One of our most popular frames, the Maxi fits a wide range of faces from smaller to medium. Soft silicone padding on the nose bridge and temple bars make this frame extremely comfortable while maintaining a sure grip. It is also available in a larger size. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Silver; Red; Black; Blue**Order No:** #RGMX30 Medium, #RGMX55 Large

Cool

Description: Sharp corners and stylish angles make the 703 a frame you can be happy about wearing. The comfortable wraparound frame comes with perforated rubber nosepieces for increased breathability and grip—an essential feature on leaded eyewear. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Grey; Black; Brown Fade

Order No: #RG-703

**Spider**

Description: The Model 230 offers women a chance to get away from standard safety frames and move towards a chic look without sacrificing safety or functionality. The 230 has a head-hugging temple design and excellent splash and side protection and is made out of a durable TR-90 nylon plastic for high impact safety and chemical resistance. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Red; Black

Order No: #RG-230

**Hollywood**

Description: The W200 is a premium protective women's eyewear frame centered on protection, performance, and a little bit of fun. Excellent splash and side protection, head-hugging temple bars with star decals, Jackie-O sunglasses-style big lenses, and a durable TR-90 frame material for impact and chemical resistance all come together to make a great pair of glasses that you can feel safe wearing and happy being seen in. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Blue; Black

Order No: #RG-W200

**Flame**

Description: The 99 is an alluring wraparound frame with a slick embossed flame decal on the temple. Lightweight, durable, comfortable, and attractive, the 99 gives you all that you could ask for and more in a safety frame. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Red; Black; Silver

Order No: #RG-99

**Speedo**

Description: The 206 utilizes a classic safety style that is solid, sturdy, and comfortable. This somewhat narrow-fitting frame is appropriate for smaller to medium faces. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Yellow/Black; Orange/Brown

Order No: #RG-206

**Wiley Romer II**

Description: The Romer II Advanced offers great coverage in a stylish frame and comes in a Matte Black finish. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Matte Black

Order No: #RG-1006F

**Wiley P17**

Description: The P-17 has been one of the most popular Wiley X frames for years. Available in Gloss Black (P-17F) and Brown Tortoise (P-17NPF), the P-17 offers ample visibility and protection in a wraparound lens style. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Tortoise; Black

Order No: #RG-P17



*Wiley X frames are available in many other styles. We show two of the most popular.
Call us about any other Wiley X style you may be interested in.*

Geek

Description: The 66 is a great standard for a multitude of medical and scientific applications. The temple bars are fully adjustable, ratcheting up and down and telescoping for a perfect fit. The large lens design ensures maximum safety and visibility. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Black; Silver

Order No: #RG-66



Astro

Description: The 250 is a standard economy safety frame which offers excellent protection, large viewing area, adjustable temple bar and built-in side shields for added protection. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Blue/Clear

Order No: #RG-250



Classic

Description: The 52 tortoise wrap is a classic in radiation protection frames, and ours now features a silicone no-slip unifi bridge for increased comfort. This frame has excellent peripheral vision with leaded glass side shields. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Tortoise

Order No: #RG-52



Goggle Wrap

Description: The RK2 is a vented goggle with neoprene foam on the inside to protect against splashing and sweat on your face. The flame-resistant elastic adjustable head band makes this frame extremely comfortable, and the large lens pocket allows for maximum viewing area. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Black

Order No: #RG-RK2



Goggle

Description: The 55 is a silicone goggle with soft foam on the inside to protect against splashing and sweat on your face. The elastic adjustable head band makes this frame extremely comfortable. The large lens pocket allows for maximum viewing area. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Black

Order No: #RG-55



Phenom

Description: You will look great wearing the Model 850 designer safety glasses. Made of flex nickel and hand made Nylon materials, this frame is lightweight, flexible and comfortable. Mono block and spring hinges add to the durability of this frame. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Pewter

Order No: #RG-850



Phenom Petite

Description: The 400 has a durable, stainless steel frame with textured grips on the temple bars to help reduce the chances of slipping off your head. This frame has a universal stylish look, spring hinges, and a durable construction. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Gunmetal

Order No: #RG-400

**Metalite Pewter**

Description: Stylish and durable, the 700 frame combines function with fashion. Made of 100% Titanium (Nickel-free to eliminate allergic reactions), the 700 comes with spring hinges for a comfortable fit. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Pewter

Order No: #RG-700

**Metalite Bronze**

Description: If you're looking for comfort and durability in a safety frame without sacrificing looks, then the 553 is for you. This sturdy frame is bronze with black accents around the rim and bridge. Couple that with the dark burgundy leopard-print on the earpieces and you've got a dark, stylish look to compliment any face. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Bronze

Order No: #RG-553

**Metalite Standard**

Description: The 554 is a durable safety frame that looks great and gets the job done. It features a classic square rim look and mildly adjustable rubber nose pads for comfort. Black and burgundy leopard-print earpieces compliment this frame for a finished combination of style, safety, and clear vision. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Black; Gold

Order No: #RG-554

**Fitover**

Description: The Fitover is a comfortable, lightweight, and durable black rubberized nylon frame that is designed to fit over prescription eyewear. These glasses will fit over almost any style of glasses with comfort. Featuring our standard 0.75mm lead equivalency Schott SF-6 HT radiation safety glass lenses.

Color: Black; Tortoise

Order No: #RG-33

