




medical products

1-800-833-9633 dexmed.com

Acrylic Arm Support



ITEM #800-0084 Small Guard, 6"H x 10"L

ITEM #91-634 Medium Guard, 8"H x 16"L

ITEM #800-0082 Large Guard, 10.5"H x 16"L

FEATURED NEW PRODUCTS

MRI COMFORT PILLOW

CFI Medical Solutions introduces the MRI Comfort Pillow. Made of high density foam and enclosed in a tough, cleanable nylon cover, the MRI Comfort Pillow is designed to comfortably support the patient's head and shoulders during procedures where the arms must be held over-head.

The MRI Comfort Pillow integrates two reinforced hand straps for the patient to hold during the procedure to increase comfort and providing a feeling of security.

Four convenient supply pouches are provided on the back of the pillow to allow for quick, easy access to cotton balls, alcohol wipes, needles, and other supplies.



Part #	Description of Part
--------	---------------------

463	MRI Comfort Pillow 8.25" x 25" x 16"
-----	--

CARBON FIBER ARMBOARDS

CFI's Carbon Fiber Armboards have been custom engineered to achieve a high degree of strength, durability and radiolucency. Two unique armboard designs are available: The "Rail Mount" armboard is designed to fit any standard accessory rail, and can also be combined with an optional accessory rail mounting unit, allowing for attachment to tables without a built in rail. The unique "Quick Release" Swivel Unit allows for easy application & removal, while providing 180° of lateral adjustment. The "Shoulder Mount" armboard includes a hexagonal base that is placed under the table pad and is held in place by the weight of the patient, allowing for quick and easy installation and removal. CFI also offers a wide range of armboard pads to fit your exact needs. By utilizing high grade materials including Conductive and Pressure Relief Covers, Memory Foam, Comfort Foam, Ultra Foam and Poly Foam, CFI is sure to have an Armboard Pad that provides the right combination of comfort and support.



Part #	Description of Part
--------	---------------------


250	Carbon Fiber Armboard Rail Mount
251	Armboard Mounting Unit
252	Armboard Pad Memory Foam with Scancoat, Fits 250
253	Armboard Pad Ultra Foam with Scancoat, Fits 250
254	Carbon Fiber Armboard Shoulder Mount
255	Armboard Pad Comfort Foam with Pressure Relief Cover, Fits 254
256	Armboard Pad Medium Density Foam with Pressure Relief Cover, Fits 254



Rail Mount





Shoulder Mount



MEDICAL IMAGING

TORSO/BODY POSITIONERS

Part #	Description of Part	Application	Alternatives	Modality
100	Torso Block - 27° 7"H x 16"W x 12"L	Used for dorsal and lumbar spine support.	-UF, -8, -20, -SCW, -SCB	GP, CT, MR, NM, US
107	Angular Bolster 7.25"H x 12.5"W x 21"L	Used for torso and aortogram procedures.	-UF, -8, -20, -SCW, -SCB	GP, CT, MR, NM, US
109	Angular Bolster 7.25"H x 12.5"W x 15"L	Specifically designed for mammograph procedures.	-UF, -8, -20, -SCW, -SCB	GP, CT, MR, NM, US
115	Bolster 5"H x 7"W x 15"L	Used to support knees and back.	-UF, -8, -20, -SCW, -SCB	GP, CT, MR, NM, US
124	Bolster 6.6"H x 9"W x 18"L	Used specifically to stabilize patient for myelogram needle insertion.	-UF, -8, -20, -SCW, -SCB	GP, CT, MR, NM, US
145	Spinal Wedge 7"H x 7"W x 21.25"L	Supports the oblique of the spine.	-UF, -8, -20, -SCW, -SCB	GP, CT, MR, NM
146	Cylinder 6"D x 12"L	Used for abdominal compression during kidney biopsies.	-UF, -8, -20, -SCW, -SCB	GP, CT, MR, NM
153	Cylinder 4.25"D x 12"L	Used for abdominal compression during kidney biopsies.	-UF, -8, -20, -SCW, -SCB	GP, CT, MR, NM



alternatives key: -UF ultra foam, -8 standard vinyl cover,
-20 heavy duty vinyl cover, -SCW scancoat® white, -SCB scancoat® black

modality key: GP general positioning,
CT CT scan, MR MRI, NM nuclear, US ultrasound



MEDICAL IMAGING GENERAL POSITIONING

All rectangle and wedge blocks are available in the following alternatives: -UF, -8, -20, -SCW, -SCB. (see page 4)
For custom wedges and rectangles, please fill out the form on page 41 and fax to 810.750.5310 for a quote.

RECTANGLE AND WEDGE BLOCKS

103 Rectangle
3"H x 8"W x 12"L



114 Rectangle
2"H x 7.5"W x 9.5"L



126 Set of 4 Rectangles
(1) .5" x 4" x 6"
(1) 1" x 4" x 6"
(1) 1.5" x 4" x 6"
(1) 2" x 4" x 6"



131 Rectangle
4"H x 10"W x 10"L



142 Rectangle
3"H x 4.75"W x 9.5"L



143 Rectangle
4.75"H x 4.75"W x 9.5"L



151 Rectangle
4"H x 18"W x 24"L



152 Decubitis Block
3"H x 14"W x 18"L



154 Set of 2, Thin Mats
.25"H x 18"W x 24"L
(only available in plain foam and -8)



169 Decubitis Block
3"H x 7.5"W x 10"L



105 Set of 2, 18.75° Wedges
3.75"H x 9"W x 8"L



106 45° Wedge
6.6"H x 6.6"W x 12"L



108 45° Wedge
5"H x 5"W x 15"L



110 Set of 2, 45° Wedges
3"H x 3"W x 7.25"L



111 45° Wedge
6.75"H x 6.75"W x 7.25"L



113 Set of 2, 12° Wedges
2.25"H x 9.5"W x 7.25"L



116 Set of 2, 34° Wedges
3"H x 4"W x 7"L



117 Set of 2, 16° Wedges
2.5"H x 7"W x 10"L



118 45° Wedge
7"H x 7"W x 7"L



120 18° Wedge
2.5"H x 8"W x 7"L



122 Set of 2, 22° Wedges
3.5"H x 7"W x 8"L



123 33° Wedge
6.6"H x 9"W x 12"L



128 13.5° Wedge
3"H x 10"W x 10"L



129 21° Wedge
4.6"H x 10"W x 10"L



130 45° Wedge
10"H x 10"W x 10"L



136 12.5° Wedge
2.25"H x 8"W x 10"L



138 21° Wedge
3.75"H x 8"W x 10"L



139 45° Wedge
7"H x 7"W x 10"L



140 27° Wedge
6"H x 10"W x 10"L



141 45° Spinal Wedge
9"H x 9"W x 30"L



144 45° Wedge
7"H x 7"W x 14"L



145 45° Spinal Wedge
7"H x 7"W x 21.25"L



158 Set of 2, 15° Wedges
3"H x 11"W x 7"L



159 Set of 2, 20° Wedges
4"H x 11"W x 7"L



191 45° Wedge
7"H x 7"W x 24"L



201 55° Wedge
11"H x 7"W x 30"L



alternatives key: -UF ultra foam, -8 standard vinyl cover, -20 heavy duty vinyl cover, -SCW scancoat® white, -SCB scancoat® black

modality key: GP general positioning, CT CT scan, MR MRI, NM nuclear, US ultrasound